

VITA ADIVA® LUTING SOLUTIONS

The material-specific luting solution – intuitive and process-proof!



VITA shade determination

VITA shade communication

VITA shade reproduction

VITA shade control

Date of issue: 07.19

VITA – perfect match.


VITA

Product description

VITA ADIVA LUTING SOLUTIONS are luting agents perfectly matched with VITA materials for bonding restorations made of feldspar, hybrid and glass ceramics, zirconia, composite and metal. The integrated material-specific coupling agents ensure dependable adhesive bonding. Clearly and easily structured, the sophisticated system supports you during

process-proof and efficient bonding, thanks to its intuitive and simple handling. All components are easy to access for you and your team, thanks to the handy set in a standard tray. VITA ADIVA LUTING SOLUTIONS are also suitable for restorative materials of other manufacturers.

ALL SOLUTIONS AT A GLANCE



VITA ADIVA® FULL-ADHESIVE

- The full-adhesive luting system with perfectly matched components for conditioning tooth substance and restorative material.
- The dual-curing, full-adhesive VITA ADIVA F-CEM luting composite is available in four shades (A2 Universal, A3, White opaque, Translucent).

VITA ADIVA® SELF-ADHESIVE*

- The self-adhesive luting system with perfectly matched components for conditioning the restorative material.
- The dual-curing, self-adhesive VITA ADIVA S-CEM luting composite is available in two shades (A2 Universal, Translucent).

VITA ADIVA® IMPLANT

- VITA ADIVA IA-CEM is a dual-curing, ultra-opaque luting composite for bonding crowns on titanium bases and individual titanium abutments.

VITA ADIVA® TEMPORARY

- VITA ADIVA TE-CEM is a dual-curing composite cement for temporary bonding. It is available in the Translucent shade.

*) VITA ADIVA S-CEM is not certified by FDA.

Benefits



Reliable

- Bonding numerous restorative materials with material-specific, matched coupling agents in a dependable way.

Process-proof

- Process-proof bonding, thanks to a clearly structured luting set with systematically arranged and color-coded components.

Efficient

- All components are easy to access, thanks to the handy, two-piece set in a standard tray.

VITA ADIVA® FULL-ADHESIVE LUTING SET



VITA ADIVA® FULL-ADHESIVE LUTING SET, Prod. No. FAFAS

| Designation | Quantity | Available as | Prod. No.* |
|--|----------|--------------------------------------|------------|
| VITA ADIVA TOOTH-ETCH** | 3 ml | Syringe + 15 application tips, blue | – |
| VITA ADIVA T-BOND I + II | 5 ml | Bottle, 2x | FATBS |
| VITA ADIVA CERA-ETCH** | 3 ml | Syringe + 20 application tips, red | FACE3 |
| VITA ADIVA C-PRIME | 5 ml | Bottle | FACP5 |
| VITA ADIVA ZR-PRIME | 5 ml | Bottle | FAZRP5 |
| VITA ADIVA M-PRIME | 5 ml | Bottle | FAMP5 |
| VITA ADIVA F-CEM, shade: A2 Universal | 5 ml | Syringe + 15 T-Mixers 1:1, brown | FAFCEMA25 |
| VITA ADIVA OXY-PREVENT | 3 ml | Syringe + 10 application tips, green | FAOP3 |
| VITA ADIVA ACCESSORIES: VITA ADIVA Microbrush, VITA ADIVA mixing trays | | | |

VITA ADIVA® FULL-ADHESIVE LUTING VERSIONS

| Designation | Quantity | Available as | Prod. No.* |
|---------------------------------------|----------|----------------------------------|------------|
| VITA ADIVA F-CEM, shade: A2 Universal | 5 ml | Syringe + 15 T-Mixers 1:1, brown | FAFCEMA25 |
| VITA ADIVA F-CEM, shade: A3 | 5 ml | Syringe + 15 T-Mixers 1:1, brown | FAFCEMA35 |
| VITA ADIVA F-CEM, shade: white opaque | 5 ml | Syringe + 15 T-Mixers 1:1, brown | FAFCEMW05 |
| VITA ADIVA F-CEM, shade: translucent | 5 ml | Syringe + 15 T-Mixers 1:1, brown | FAFCEMT5 |


VITA ADIVA® FULL-ADHESIVE ACCESSORIES

| Designation | Quantity | Prod. No.* |
|-------------------------|----------|------------|
| VITA ADIVA Microbrush | 25 pcs. | FAMB25 |
| VITA ADIVA mixing trays | 25 pcs. | FAAT25 |
| VITA ADIVA T-Mixer 1:1 | 15 pcs. | FATM11 |

*) All VITA ADIVA components can be (re)ordered individually using the Prod. No. specified.

**) Note: Refill delivery form for VITA ADIVA TOOTH-ETCH: 2 x 3 ml syringe, Prod. No. FATE3; VITA ADIVA CERA-ETCH also available as 6 ml bottle, Prod. No. FACE6

VITA ADIVA® SELF-ADHESIVE LUTING SET



VITA ADIVA® SELF-ADHESIVE LUTING SET, Prod. No. FASAS

| Designation | Quantity | Available as | Prod. No.* |
|--|----------|--------------------------------------|------------|
| VITA ADIVA CERA-ETCH** | 3 ml | Syringe + 20 application tips, red | FACE3 |
| VITA ADIVA C-PRIME | 5 ml | Bottle | FACP5 |
| VITA ADIVA ZR-PRIME | 5 ml | Bottle | FAZRP5 |
| VITA ADIVA M-PRIME | 5 ml | Bottle | FAMP5 |
| VITA ADIVA S-CEM, shade: A2 Universal | 5 ml | Syringe + 15 T-Mixers 4:1, brown | FASCEMA25 |
| VITA ADIVA OXY-PREVENT | 3 ml | Syringe + 10 application tips, green | FAOP3 |
| VITA ADIVA ACCESSORIES: VITA ADIVA Microbrush, VITA ADIVA mixing trays | | | |

VITA ADIVA® SELF-ADHESIVE LUTING VERSIONS

| Designation | Quantity | Available as | Prod. No.* |
|---------------------------------------|----------|----------------------------------|------------|
| VITA ADIVA S-CEM, shade: A2 Universal | 5 ml | Syringe + 15 T-Mixers 4:1, brown | FASCEMA25 |
| VITA ADIVA S-CEM, shade: translucent | 5 ml | Syringe + 15 T-Mixers 4:1, brown | FASCENT5 |

VITA ADIVA® SELF-ADHESIVE ACCESSORIES

| Designation | Quantity | Prod. No.* |
|-------------------------|----------|------------|
| VITA ADIVA Microbrush | 25 pcs. | FAMB25 |
| VITA ADIVA mixing trays | 25 pcs. | FAAT25 |
| VITA ADIVA T-Mixer 4:1 | 15 pcs. | FATM41 |

*) All VITA ADIVA components can be (re)ordered individually using the Prod. No. specified.

**) Note: VITA ADIVA CERA-ETCH is also available as 6 ml bottle, Prod. No. FACE6

VITA ADIVA® IA-CEM and TE-CEM



VITA ADIVA® IA-CEM SET Prod. No. FAIACEMSET

| Designation | Quantity | Available as | Prod. No.* |
|---|----------|----------------------------------|------------|
| VITA ADIVA IA-CEM, shade: ultra opaque | 5 ml | Syringe + 15 T-Mixers 1:1, brown | FAIACEMU05 |
| VITA ADIVA CERA-ETCH | 6 ml | Bottle | FACE6 |
| VITA ADIVA C-PRIME | 5 ml | Bottle | FACP5 |
| VITA ADIVA ZR-PRIME | 5 ml | Bottle | FAZRP5 |
| VITA ADIVA IA-CEM ACCESSORIES: VITA ADIVA T-Mixer, VITA ADIVA Microbrush, VITA ADIVA mixing trays | | | |

VITA ADIVA® IA-CEM ACCESSORIES


| Designation | Units | Prod. No.* |
|-------------------------|-----------|------------|
| VITA ADIVA T-Mixer 1:1 | 15 pieces | FATM11 |
| VITA ADIVA Microbrush | 25 pieces | FAMB25 |
| VITA ADIVA mixing trays | 25 pieces | FAAT25 |

VITA ADIVA® TE-CEM

| Designation | Quantity | Available as | Prod. No. |
|--------------------------------|----------|----------------------------------|-----------|
| VITA ADIVA TE-CEM, translucent | 5 ml | Syringe + 15 T-Mixers 4:1, brown | FATECEM5 |

*) All VITA ADIVA components can be (re)ordered individually using the Prod. No. specified.

SHADES AND TYPES









MATERIAL-SPECIFIC BONDING RECOMMENDATIONS

| VITA CAD/CAM material | Material class | Type of bond / luting-composite | | | | |
|-----------------------|---|---------------------------------|-------------------|--------------------|----------------|--------------------|
| | | adhesive | self-adhesive | adhesive | conventional | temporary |
| | | VITA ADIVA® F-CEM | VITA ADIVA® S-CEM | VITA ADIVA® IA-CEM | | VITA ADIVA® TE-CEM |
| VITABLOCS | feldspar ceramics | ● | ● ¹ | ○ | X | X |
| VITA ENAMIC | Hybrid ceramic | ● | ● ¹ | ● | X | X |
| VITA SUPRINITY PC | Zirconia-reinforced Lithium silicate ceramic | ● | ● | ● | ○ ² | X |
| VITA YZ T | Zirconia | ● | ● | ● | ● | ● |
| VITA YZ HT | Zirconia | ● | ● | ● | ● | ● |
| VITA YZ ST | Zirconia | ● | ● | ● | ● | ● |
| VITA YZ XT | Zirconia | ● | ● | ● | ● | ● |
| VITA CAD-Temp | Composite | ○ | X | ● | ○ | ● |
| - | NPM, PM, titanium | ● | ● | ● | ● | ● |

● recommended ○ possible X not possible

TABLE OF INDICATIONS

| | VITA ADIVA® F-CEM | VITA ADIVA® S-CEM | VITA ADIVA® IA-CEM | VITA ADIVA® TE-CEM |
|---|--------------------------------|--------------------------------|--------------------------------|---------------------|
| | Full-adhesive luting composite | Self-adhesive luting composite | Full-adhesive luting composite | Temporary composite |
|  | ● | ● | ● | ● |
|  | ● | X | X | X |
|  | ● | X | X | X |
|  | ● | ● | ● | ● |
|  | ● | ● | ● | ● |
|  | ● | ● | ● | ● |


● recommended X not possible

Please note:

¹) In the case of silicate ceramics and hybrid ceramics, self-adhesive bonding composites may only be used for crowns.

²) Conventional cementation materials (Ketac CEM, 3M ESPE; Vivaglass CEM, Ivoclar Vivadent, GC Fuji I, GC Dental) may be used exclusively for crowns on prepared natural teeth. The preparation requires retentive surfaces for conventional cementation, and - based on the preparation guidelines - the anatomical shape must be reduced while adhering to the minimum layer thickness indicated.

More information about VITA ADIVA LUTING SOLUTIONS is available at www.vita-zahnfabrik.com/adiva



Please note: Our products must be used in accordance with the instructions for use. We accept no liability for any damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product. Date of issue of this information: 07.19

After the publication of this information for use any previous versions become obsolete. The current version can be found at www.vita-zahnfabrik.com

VITA Zahnfabrik has been certified and the following products bear the CE mark:
CE 0482

VITA ADIVA® F-CEM, VITA ADIVA® S-CEM, VITA ADIVA® IA-CEM, VITA ADIVA® TE-CEM, VITA ADIVA® TOOTH-ETCH, VITA ADIVA® T-BOND I + II, VITA ADIVA® C-PRIME, VITA ADIVA® ZR-PRIME, VITA ADIVA® M-PRIME

Harvard Dental International GmbH has been certified in accordance with the Medical Device Directive and the following products bear the CE mark:
CE

VITA ADIVA® Microbrush®, VITA ADIVA® T-Mixer™

VITA Zahnfabrik has been certified in accordance with the Medical Device Directive and the following product bears the CE mark:
CE

VITA ADIVA® OXY-PREVENT

Microbrush® is a trademark of MICROBRUSH® INTERNATIONAL.
T-Mixer™ is a trademark of Sulzer Mixpac AG.

The products/systems of other manufacturers mentioned in this document are registered trademarks of the respective manufacturers.

VITA

VITA Zahnfabrik H. Rauter GmbH & Co. KG
Spitalgasse 3 · D-79713 Bad Säckingen · Germany
Tel. +49(0)7761/562-0 · Fax +49(0)7761/562-299
Hotline: Tel. +49(0)7761/562-222 · Fax +49(0)7761/562-446
www.vita-zahnfabrik.com · info@vita-zahnfabrik.com
 facebook.com/vita.zahnfabrik